

# An Assistive Wheelchair-Mounted Frisbee Launcher



WPI


David Barsoum, Alexander Kaneko, Alex Chen, Nathan Lam

## Problem

Due to mobility issues in the shoulder and elbow, those who use wheelchairs cannot throw a frisbee autonomously to play with their dog, friends, or family.

## Features

- Launches frisbees up to 55 feet
- Supports a large range of frisbee sizes
- Switch adapted and accessible
- Tested with precision machining
- Made with easily obtainable materials (wood, cardboard, acrylic, metal)



Final Prototype

## Objective


Design a wheelchair-mounted projectile-launching device to assist clients with throwing frisbees.

## Design Studies

- 1 Clip Testing**  
Measured launching distances of prototype with and without clip to improve stability without sacrificing power.
- 2 Spring Testing**  
Observed the impact of spring constant on launch distance by testing the launcher with springs of varying levels of strength.
- 3 Mount Stability Study**  
Tested the stability of two mounting systems by considering how well they supported weight.


## Designs

1




Initial cardboard launcher

2




Manually-operated spring-powered launcher

3



Motorized launcher with plastic gears

4



Motorized launcher with aluminum gears



SOLIDWORKS

# Operation of Device

1

Mount the launching device to the wheelchair.

2

Turn the electronic control system switch to "on"

3


Insert the frisbee into the cardboard arm.

4

Make sure that the direction of the motor is correct:  
"forward" for launching and "reverse" to un-tension

5

Press and hold the "LAUNCH" button until the arm is flung  
and the frisbee is launched forward.



## Safety




Wear safety glasses when launching




Extend



Extend the launcher away from you  
using the slides when launching.



Un-tensioned Spring



Tensioned Spring

When not in use, the launcher must be  
stored in an un-tensioned state. Turn the  
electronic control system off and  
disconnect the battery.



SOLIDWORKS